



MICRON[®]

Computer Requirements

On the web: <https://phoenixmicron.com> **Email us at:** sales@phoenixmicron.com

MICRON[®] SE and **MICRON Review Station** each allow you to supply your own computer system. The chart below will help ensure that you choose the appropriate computer for your application.

Optionally you can purchase a preconfigured MICRON Control Unit directly from Phoenix-Micron or your distributor.

Requirements	MICRON [®] SE Imaging System	MICRON [®] Review Station
Processor (CPU)	Intel [®] Core™ i7 (13 th Gen +)	Intel [®] Core™ i7 (13 th Gen +)
RAM	16 GB DDR4	16 GB DDR4
Graphics Card (GPU)	Intel [®] UHD Graphics 730	Intel [®] UHD Graphics 730
SSD	512GB NVMe PCIe Gen 3	512GB NVMe PCIe Gen 3
HDD (for image storage)	2TB	2TB
Operating System	Windows 10 Pro (64-bit) Windows 11 Pro (64-bit)	Windows 10 Pro (64-bit) Windows 11 Pro (64-bit)
Additional Hardware Requirements		
USB Ports	Min. 3x USB 3.2	Optional
Ethernet Port	1x Gigabit	Optional
Display Port	Display port or HDMI (monitor dependent)	Optional

